



Software Release Note and Revision History Example

The training document example is fictitious and should only be applied in the appropriate context with oversight by regulatory and software professionals with direct knowledge and experience with software configuration management. The document should not be used as a cookbook or taken literally without knowledgeable evaluation of current interpretations and enforcement.

While SoftwareCPR® attempts to ensure the accuracy of information presented, no guarantees are made since regulatory interpretations and enforcement practices are constantly changing and are not entirely uniform in their application.

Disclaimer of Warranties: The information is provided AS IS, without warranties of any kind. CPRLLC does not represent or warrant that any information or data provided herein is suitable for a particular purpose. CPRLLC hereby disclaims and negates any and all warranties, whether express or implied, relating to such information and data, including the warranties of merchantability and fitness for a particular purpose.

For additional support, contact SoftwareCPR at office@softwarecpr.com or call +1 781-721-2921

Acme Device System Version History

The Acme Device System software components are described in the ADD-10-002R01 Architectural Design Description.

Current Production Release:

Date	System Revision	Description of Change												
05 Mar 2019	2.1.0	New features and fixes: <ul style="list-style-type: none"> • Added average pump rate trend display. • Bug with inability to save the changes to the arcwires has been fixed. • Opening more than one patient case within the same DIBS instance is no longer supported. • Corrected open anomalies: <ul style="list-style-type: none"> ○ DF-0089 – Corrected issue with host comm status bit reflecting disconnect. 												
		<table border="1"> <thead> <tr> <th>Custom Components</th> <th>Revision</th> </tr> </thead> <tbody> <tr> <td>Blood Pump</td> <td>1.1.0</td> </tr> <tr> <td>Venous Line Clamp</td> <td>1.0.5</td> </tr> <tr> <td>Air Detector</td> <td>1.1.2</td> </tr> <tr> <td>UI</td> <td>2.1.0</td> </tr> <tr> <td>Host Comm</td> <td>1.5.1</td> </tr> </tbody> </table>	Custom Components	Revision	Blood Pump	1.1.0	Venous Line Clamp	1.0.5	Air Detector	1.1.2	UI	2.1.0	Host Comm	1.5.1
		Custom Components	Revision											
		Blood Pump	1.1.0											
		Venous Line Clamp	1.0.5											
		Air Detector	1.1.2											
		UI	2.1.0											
		Host Comm	1.5.1											
		SOUP/OTS Software Components												
		Xeno Library	6.3.0											
Qt	5.12 LTS													
FastMath Library	1.55													

Revision History:

Date	System Revision	Description of Change												
22 Nov 2018	2.0.0	New features and fixes: <ul style="list-style-type: none"> • Added new graphical user interface for overall system control. • Added smooth start feature for initial pump start. • Corrected open anomalies: <ul style="list-style-type: none"> ○ DF-0101 – ... ○ DF-0099 – ... ○ DF-0091 – ... 												
		<table border="1"> <thead> <tr> <th>Custom Components</th> <th>Revision</th> </tr> </thead> <tbody> <tr> <td>Blood Pump</td> <td>1.1.0</td> </tr> <tr> <td>Venous Line Clamp</td> <td>1.0.5</td> </tr> <tr> <td>Air Detector</td> <td>1.1.2</td> </tr> <tr> <td>UI</td> <td>2.0.18</td> </tr> <tr> <td>Host Comm</td> <td>1.5.0</td> </tr> </tbody> </table>	Custom Components	Revision	Blood Pump	1.1.0	Venous Line Clamp	1.0.5	Air Detector	1.1.2	UI	2.0.18	Host Comm	1.5.0
		Custom Components	Revision											
		Blood Pump	1.1.0											
		Venous Line Clamp	1.0.5											
		Air Detector	1.1.2											
		UI	2.0.18											
		Host Comm	1.5.0											
		SOUP/OTS Software Components												
		Xeno Library	6.1.0											
Qt	5.11.3													
FastMath Library	1.55													

...

...

< Generally one would expect to see list any version released to the field or used for verification and validation activities. >

For additional support, contact SoftwareCPR at office@softwarecpr.com or call +1 781-721-2921