


spectral evolution

ILM-550 Tungsten-Halogen Illuminator Product Sheet

- ▶ **50 watt Tungsten**-Halogen source with reflector for bright, uniform illumination
- ▶ **Insulated handle, swivel mount and stand for easy positioning** on bench surfaces and optical tables
- ▶ **Includes 150 watt, 12 volt power supply**— ample power for higher wattage bulbs if needed!



SPECTRAL EVOLUTION'S ILM-550** Tungsten-Halogen illuminators feature user-replaceable bulbs in a standard socket layout— you can upgrade your illuminator to a 75W or 100W bulb for more power if needed. Each unit comes **with a 3" diffuser mounted in a removable screw mount for use with filters!

Contact your local SPECTRAL EVOLUTION representative today for additional information!



1 Canal Street ◊ Unit B1
Lawrence, MA 01840 USA
Tel: 978 687-1833 ◊ Fax: 978 945-0372
Email: sales@spectralevolution.com
www.spectralevolution.com